

# PROGRESSIVE DENTAL SUPPLY



[www.SmileStream.com](http://www.SmileStream.com)



SMILE & SAVE **MORE** WHEN YOU BUY INVENTORY SUPPLIES WITH THESE DISCOUNTS!



# FREE UPS SHIPPING

ON ALL ORDERS OVER

**\$500**<sup>USD</sup>

USA ORDERS  
(GROUND SHIPPING ONLY)

**\$650**<sup>USD</sup>

INTERNATIONAL  
ORDERS



## BETTER THAN ROTH BRACKET SYSTEM

**35% OFF**

50 OR MORE METAL

**15% OFF**

10 OR MORE CERAMIC

Use code: BTR METAL

Use code: BTR CERAMIC

MIX AND MATCH 50 BTR METAL COLORS OR 10 BTR CERAMIC COLORS AND RECEIVE A DISCOUNTED PRICE!

(CANNOT MIX METAL AND CERAMIC)

## TABLE OF CONTENTS

Brackets .....	1
Bands .....	3
Archwires .....	6
Adhesives.....	8
Separators.....	11
Tie In Materials.....	12
Chain & Coil.....	13
Hooks, Cleats & Expansion Products.....	15
Instruments.....	17
Headgear.....	25
Elastics.....	27
Peripheral Aligners.....	28
Waxes.....	30
Patient Mgmt & Records...31	
Packaging Templates.....33	
Bracket Wheels.....36	
Band Fitting Kits.....38	
Skeletal Anchorage.....38	
Education & Office Kits....43	
Books, Brochures &.....45	
Office Marketing Materials	



\*Some items in catalog may not be available from the PDS Spain Warehouse or outside the United States.

## BRACKETS

The IP Appliance® system has a variety of bracket prescriptions available, so that you can be more effective in your treatment of each case. By selecting the perfect appliance for each case, your ortho will be more efficient, and has greater consistency among your patients. This is the logic behind the IP Appliance®, generally referred to as a “customized appliance”.

All IP brackets are cut with .022 slots, with the exception of a few .020 slot (SLa and SLi) brackets on the incisors. This allows the application of more rigid archwires, which reduce unwanted tooth movements (especially in large extraction spaces) and level the curve of spee.



### Bracket Variations Available

**Standard Roth:** The standard Roth type bracket is the standard in the Orthodontic industry.

**Torque (La, Li):** In class III cases, extra labial root torque should be added to prevent excess incisor proclination. The added labial root torque can help improve the inclination and esthetics of the incisors, allowing for more non-extraction treatments. Li brackets are standard in dental and skeletal deep bite cases where the incisors are retroclined. In class II cases, extra lingual (Li) root torque added to the upper incisors and extra labial root torque added to the lower incisors will result in better, more efficient treatment.

**Rotation (D, M):** Rotation Brackets provide 4 degrees of overcorrection to eliminate finishing bends and to improve the retention experience.

**Combo (MLa, MLi, DLa, DLi):** Combo brackets combine both torque and rotation prescriptions.

**“Super” Torque:** These brackets have a reduced slot size of .020.

### Other Bracket Options Available

**Rs:** Standard Roth bonded bracket with a small bonding pad. Used on teeth with short clinical crowns.

**Pe:** Partial erupted bonded bracket to be applied to second molars that are not fully erupted through the tissue. This bracket is especially needed in the upper 7-extraction diagnosis to prevent super-eruption of the lower molar.

**CIIE:** “Bonded” convertible tube with added mesial rotation.

**Tie:** Single tie-wing bracket to be used when both central incisors have a “mesial” rotation.

**Th:** Low profile bracket to be applied to deep bite cases.

**Up:** Molar uprighting bonded bracket. This convertible buccal tube has added mesial rotation, and added lingual root torque to correct lingual crown inclination.

### Bracket Kits

#### Metal Bracket Repair Kit (011-004):

The Metal Bracket Repair Kit consists of two metal brackets of every variation. The brackets are used as back-ups in the clinic in case a patient loses brackets, or if you need to alter your prescription.

#### AM Bracket Repair Kit (011-005)

The “AM” (Absolute Minimum) Metal Bracket Repair Kit consists of one metal bracket of every variation. AM Kits are designed to supply you with a small amount of inventory.





## BANDS

In the IP Appliance® band system, fitting bands are used to determine the size of the band needed for your patient.

### Fitting Bands

Fitting bands come with eyelets to assist in seating. Our fitting band sizes exactly match the size of the treatment bands that we will weld the bracket to.



### Complete Band Fitting Kit (015-010-00)

The Complete Band Fitting Kit comes with molar fitting bands for all 4 quadrants, along with the upper and lower bicuspid. The 4 molar quadrants contain 2 fitting bands of every size from 1-34 and 1 band each of upper and lower bicuspid sizes from 1-24.



### AM Band Fitting Kit (015-010-01)

The concept of the AM Band Fitting Kit is similar to that of the Complete Band Fitting Kit, except that we have removed the bicuspid bands and are given one band of each size. Kit comes with all 4 quadrants containing 1 molar fitting band of every size from 1-34.

## Treatment Bands

Treatment Bands are the bands used in the actual case. IP Appliance® Bands have options for attachments on both the lingual and labial, as well as tip welding variations for the labial attachments.

### Lingual Attachments




**C:** Lingual Cleat

**Sh:** Lingual Sheath Goshgarian type 036x072 tube with a ball hook and lock on the distal used for Transpalatal arch (TPA).



Laser etched labeling on the side of the band indicate quadrant and the band size. For this band, 19UR is the size 19 band in the Upper Right Quadrant.

IP Appliance® Bands use a simple numbering system to allow the user to easily understand the attachments on each band. The item number is as follows:



## Bicuspid Bands

Variations available in Bicuspid Bands:

- **R2:** Standard Roth type bracket
- **M2:** Mesial Rotation bracket
- **D2:** Distal Rotation bracket



Welding Variations

- **Tip D**
- **Tip M**

## Labial Attachments



**R2:** Convertible Roth type bracket welded to a band.



**D2:** Convertible buccal tube without distal offset to a band. This variation is used on the upper first molars (6s) in nonextraction class III treatment.



**Up2:** Molar uprighting bracket welded to a band. This convertible buccal tube has added mesial rotation, a custom weld for distal crown tipback, and added lingual root torque to correct lingual crown inclination.



**R2 (Second Molar):** For the standard Roth type bracket on the second molar, the buccal tube is not convertible.

**R2HG:** Headgear tube attached to a Roth bracket to allow for the use of headgear facebow.

**Lip2:** Lip bumper tube with a standard Roth archwire slot. This is the variation to use when a lip bumper is applied for lip habits, maintaining E space to gain archlength, or for added lower arch anchorage.

**Pig2:** This variation has a double archwire tube to allow the use of two archwires in the upper arch at the same time, plus a third headgear tube.

**Li (Second Molar):** This variation is for the lower 2nd Molar. It has added lingual root torque to be used to upright teeth with the crown “dumped” to the lingual.



## Labial Welding Variations

**CIIE2:** Convertible Buccal tube custom welded to a band with added (distal crown) tipback, added mesial rotation.

**R2TipD:** This is a welding variation where the standard Roth convertible buccal tube is welded to the band, so that the crown will tip 5 degrees distal when a straight wire is engaged.

**R2TipM:** This is a welding variation of the standard Roth convertible buccal tube where the crown will tip mesial 5 degrees when a straight wire is engaged.



## ARCHWIRES

### Start Wires

**012 NiTi (030-0079):** Nickel Titanium wire used for level and alignment. The flexibility also makes it easy to tie into severely crowded arches. The low force is ideal for the early stages of treatment. Sold in packs of 10.

**014 NiTi (030-0081):** 014 Niti can be used interchangeably with 012 Niti. 014 is a slightly stronger wire. Sold in packs of 10.



### Treatment Wires

Shapes Available: **Ovoid** (Small, Medium, Non-Extraction #1, Non-Extraction #2), **Tapered** (Small, Medium, Non-Extraction), **Square** (Medium, Large)

Sizes Available: **Stainless Steel** (.014, .016, .018, .020, .019x.025, .021x.025) and **NiTi** (.016, .018, .018x.025 heat activated, .019x.025)



### Looped Archwires

In the IP System, Progressive also offers looped archwires in Keyhole Looped and TLoop variations.

**Keyhole Looped:** 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 42mm, 44mm

**T-Loop:** 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 42mm, 44mm, 46mm



### Archwire Repair Kit (030-000)

The kit comes with two (2) replacement archwires of every size and shape offered, packaged in four easy to use binders (one for uppers, and one for lowers), one binder for both the Keyhole and T-loop Wires. **Binders also available individually.**





## ADHESIVES



### Reverse Curve Wire (030-0080)

The reverse curve wire, is used for cases with deep bite and deep curve of spee. Sold in packs of 10.



### .016x.016 SS Archwires (030-021,022)

These square, stainless steel archwires are used to fabricate utility arches. Sold in packs of 10.

**Available for Upper & Lower.**



### NiTi Utility Arch (030-055,060,065)

The anterior segment is made from nickel-titanium wire, allowing for alignment at the same time as intrusion and/or advancement from the utility arch. Sold in packs of 3.

**Offered in three sizes: 28mm, 34mm, 38mm.**



### .036 SS Straight Wire (030-0060)

Used for custom fabrication of palatal bars (TPA), usually in a straight form. Sold in packs of 10.



### Arch Marking Pencil (030-036)

For intraoral marking of the archwire where bends need to be positioned. Essential for positioning finishing bends, loop position or placement of steps. Sold in packs of six pencils.

### Bond A Braid Lingual Retainer Wire (030-075)



This is a dead soft 8-braided wire for permanent retention and is easy to adapt. Sold in packs of 8.

### Assure™ Universal Bonding Resin (6cc) (019-013)



Assure Universal Bonding Resin will increase adhesion to normal, atypical, dry or slightly contaminated teeth. This universal bonding agent bonds to metal restorations without additional primers. No drying or light curing required and can be used with any light or chemical cure system.

Resin will bond to the following surfaces without additional primers: fluorosed, deciduous, bleached or wet enamel, dentin, metal surfaces and composite restorations. **Not available in PDS Spain.**

### Bracket-Lock (019-020)

This is the light cure adhesive used as the standard for repositioning or rebonding brackets, the initial bonding of molar brackets placed by direct or indirect techniques. Bond strength is at the maximum level immediately upon light curing.

**Bracket lock is sold as kit or individual components: Resin (019-021), Adhesive Syringe (019-022), and Pro-Etch 38% (019-023).**



### Transbond XT (019-070)

Transbond XT Light Cure Adhesive bonds metal and ceramic brackets to tooth surfaces. This adhesive uses light cure adhesive technology to provide you with additional working time to ensure accurate bracket placement. Sold with 2 syringes.

**Transbond XT Bracket Primer (6ML) (019-071)**

**Transbond XT Bracket Adhesive Only (4 syringes) (019-072)**





## Pro-Lock Band Cement Kit (020-021) Economy Kit (020-023)

This is a dual cure adhesive with a slow self-set time to allow for many bands to be seated at once, and light cure to finish the set on command.



## Excel Adhesive Kit (Herbst Appliance) (020-007)

For use in the cementation of bonded Rapid Palatal Expanders (RPE).




## Transbond Plus Band Cement Kit (020-013)

Transbond Plus is a single paste band cement. The blue color allows you to see the excess cement for easy removal, before you cure.

This kit consists of 5 paste syringes (4gr each), 50 syringe tips and instructions.



## Auxiliares



## Mixing Pads (019-036)

The Mixing Pad aids in avoiding cross contamination of materials.



## “Etch Rite” Etching Gel (020-003)

Etching Gel used as an alternative to liquid etch. Gel has the advantage of staying in place, is color coded in blue, and does not drip. Sold in a 1.2 ml. Syringe (38% Phosphoric Acid). **Not available in PDS Spain.**



## Porcelain Conditioner (020-004)

Porcelain Primer is a silane coupling agent used to improve bonding between porcelain crowns and the metal bracket.



## Application Brushes (020-009)

Brushes are used for etch liquid and activator liquid application. 2" in length. Sold in packs of 25.



## Spatulas (020-010)

Sold in packs of 100.



### Flexi Bonding Retractors (020-025,026)

This disposable product has comfortable velcro straps, which allow for patient comfort as the operator gains direct bonding access to the buccal surfaces of the bicuspids, and in some cases, the first molars.



In most cases, both upper and lower arches may be bonded at the same time when using this product.

These retractors are also useful for maintaining isolation during indirect bonding, especially when the molars are to be bonded.



### Diamond Separating Disk (050-007)

Diamond disks (Double-sided) are used to create space between teeth, typically when a separator is lost or there is not enough space for two bands. Disks are .010mm. **Not available in PDS Spain.**



### Neet Springs (050-010,015)

Nickel-titanium separating springs are best suited for separation of teeth in adult patients and overtight contacts. Less force is exerted on the teeth, resulting in significantly less pain than elastomeric separators.

The springs are available in "Molar" and "Bi-Molar" sizes. Sold in packs of 25.

## SEPARATORS



### Elastic Separators - Blue (050-001)

Latex Free. Elastomeric separators are impregnated with radio-opaque material to make them visible on x-rays. Sold in packs of 1000.



### 020-025 Brass Separating Loops (050-005,006)


Brass wire is used to maintain space already created by another type of separator while waiting for the arrival of bands from the IP Appliance or the lab. Sold in packs of 100.

## TIE IN MATERIALS



### Steel Ligatures (040-001)

010 Steel ligatures are used when a positive seat is needed to deflect the archwire into the bracket slot. Sold in packs of 1000.



### Shorty Lig Twists (040-002)

The Shorty Ligature twists are preferred by some operators since they may be placed and tightened with only one hand. The shorty twist wire is not quite as strong as the preformed .010 steel ligatures. Sold in packs of 500.



### Colored Alastic Ligatures

Kids and adults love the 26 different colors to choose from. Sold as 24 ligatures per stick, 40 sticks per bag.

Assorted Color Package (040-032)

Individual Colors (040-022 - 040-070)

## CHAIN & COIL



### Elastic Chain (040-009, 010, 012, 013)



#### Clear Chain II (040-009)

Sold in 15' spool. Stocked in spools with the standard material having one millimeter of material between the modules.



#### Gray Chain II (040-010)

Sold in 15' spool. Gray chain II is applied in cases to maintain space closure or where small amounts of space need to be closed.



#### Gray "Super" Chain I (040-012)

Sold in 15' spool. This material can be used wherever space needs to be closed.



#### Clear "Super" Chain I (040-013)

Sold in 15' spool. Clear Chain and Gray Chain exhibit similar force properties, except that degradation of force occurs even more rapidly.



### Open Coil (.030 Stainless Steel)

(040-016)

3ft Spool. Open coil (.030) has a typical "spring" look and is used to open space between teeth.



### Closed Coil (.030 Stainless Steel)


(040-017)

3ft Spool. Closed coil (.030) is tightly woven coil used to maintain space between brackets.

### Nickel-Titanium Open Coil

(040-018,018A)

Nickel Titanium Open Coil is very easy to manipulate and has ideal force properties for moving one tooth.



(040-018) - 15mm (Packs of 5)

(040-018A) - 7" (Packs of 3)

### Nickel-Titanium Closed Coil

(040-020, 021)


This product is ideal for cuspid retraction and closing extraction space.



9mm Niti Closed Coil (040-020) and 12mm Niti Closed Coil (040-021) sold in packs of 5.

### NiTi Force Markers (040-019)

The NiTi Force Marker is a plastic ruler with the steps premeasured to assure the correct activation for the NiTi Closed coils. One side of the ruler is for 9mm closed coils, and the other side is for 12 mm coils. Sold in packs of 10.



## HOOKS, CLEATS & EXPANSION PRODUCTS



### Kobayashi Hooks (040-015)

Kobayashi Hooks are used to provide an attachment for elastics. Sold in packs of 100.

### .014 Koby Twists (040-002A)

Short Kobayashi Twists used with the instrument Shorty Twist Holder. Sold in packs of 100.

### Crimpable Split Hooks (090-006)

Split Hooks are used to reactivate compressed open coil. Sold in packs of 20.

### Crimpable Split Stops (090-007)

Split Stops are also used for spacing purposes. Sold in packs of 20.

### Warren Torque Springs (090-008)

Warren Springs are used for torquing of individual teeth. Sold in packs of 10.

### Bondable Metal Buttons (090-015)

Bondable Metal Buttons with rounded pads are used on the lingual of bicuspid to obtain an attachment for power chain.

Bondable buttons are also used for unerupted teeth (such as impacted cuspids or partially erupted teeth) to provide force application.

Sold in packs of 10.

### Bondable Metal Cleats (090-004,005)

Bondable Metal Cleats are used on cuspids for attachment of lingual power products. The cleat design is more low-profile and interferes less with the occlusion than buttons.

Cleats may also be bonded to the linguals of incisors, as a reminder to the patient receiving myofunctional training of anterior tongue thrust habits.

Available in Small (090-004) and Large (090-005). Sold in packs of 10.

### Saif Springs (090-025)

10mm in length. Saif Springs are used for Class II or Class III applications when patient cooperation is lacking. Sold in packs of 2. **Not available in PDS Spain.**

### Fluted Finishing Bur (090-026)

Used to remove the adhesive remaining on the tooth after bracket debonding.

### Safety Expansion Key (090-080)

Used for RPE Activation. This product is the answer to patients swallowing keys!

### Cuspid Eruption Appliance (090-082)

Bondable eyelet with a chain attached used for retrieval of impacted cuspids. Allergic reactions minimized due to the chain being 14K gold.



## Goshgarian Transpalatal Bars (090-030,035,040,045,050,055,060,065)

Goshgarian Transpalatal Bars are used for:

- expansion of the upper arch
- derotating (mesial-buccal) upper first molars to gain arch length and Class II correction
- torquing of upper molars, anchorage enhancement
- intrusion of upper molars
- distalization of molars

Sold in packs of 5.

Available in the following sizes:

**090-030** - Full Set of 35 Palatal Arches (5/ea)

**090-035** - Palatal Bar Size 1- 34 mm

**090-040** - Palatal Bar Size 2- 37 mm

**090-045** - Palatal Bar Size 3- 41 mm

**090-050** - Palatal Bar Size 4- 45 mm

**090-055** - Palatal Bar Size 5- 47 mm

**090-060** - Palatal Bar Size 6- 49 mm

**090-065** - Palatal Bar Size 7- 52 mm



## INSTRUMENTS

### Band Pusher (060-001)

Band Pusher is used for seating and burnishing bands into their final position.

### Band Bite Stick (060-002)

Used to fit the final seat to bands by having patient bite down on the stick.



### Elastic Separating Pliers (060-004)

Used for placing elastomeric separators.



### Fast Elastic Positioner (060-005)

Used for quick and efficient placement of Pro-Ties (elastomeric ligatures).



### Mathieu Pliers (060-006)

Used for tightening ligature pigtails.



### Shorty Twist Holder (060-007)

This product is ideal for single operator reties. **Not available in PDS Spain.**



### Pin and Ligature Cutter (060-008)

Used for cutting ligature wire only.



### Hard Wire Cutter (060-009)

Designed to cut any size wire.



### Light Wire Pliers (060-010)

Used for general wire bending and the insertion of archwires into tubes.



### Jarabak Pliers (060-011)

Used for loop bending in small archwires.



### Distal End Cutter with Safety Hold (060-012)

Used to cut all archwires distal to the molar tube, or to cut stainless steel wires in your hand.



### Angle Howe Pliers (060-013)

The Angle Howe Plier has flat beaks that may be used to flatten the U loop of transpalatal bars as an intraoral adjustment.



### Posterior Band Remover (060-014,028)

Quick and efficient removal of cemented bands. (060-028) - Replacement Plastic Tips (3pk)



### Bracket Removing Pliers (060-015)

This tool is very efficient for the removal of both anterior and posterior bonded brackets.



### Three-Prong Facebow Pliers (060-016)

Use to contour the heavy wire used in head-gear facebows.



### Three Prong Adjustment Pliers (060-017)

Used for constricting loops on Hawley Bows, adjusting the Transpalatal Bar away from the palate, and reactivating quad helix appliances.



### Loop Forming Pliers (060-019)

Forms vertical loops and crimps at each end of an angle. Good for making loops for tieback or elastic placement.



### Cap Removal Tool (060-022,021)

This tool is necessary when archwire bends are needed between the molars or when inserting heavy rectangular archwires through successive molar tubes.

(060-021) - Tool Tips Only



### Tweed Loop Forming Pliers (060-025)

Intraoral adjustment of utility arches is made easy with the Tweed Pliers.

## Ligature Director (060-029)



Used to direct chilled heat activated nickel-titanium archwires into the bracket slots during retie procedures, it is also used to push steel ligatures into difficult areas to engage the bracket tie-wings, and to push open coil when placing split hooks to reactivate.

The other end of the ligature director has a serrated end to push pigtails under the bracket after cutting. Essential instrument for efficient tie-in of archwires.

## Weingardt Pliers (060-031)



Commonly used to make cinchback bends for loop activation or securing archwires.

## Lingual Arch Forming Pliers (060-035)



Forms doubleback bend in .036 wire for use in Goshgarian transpalatal arch sheaths. This plier is also used whenever a machine-made doubleback bend in preformed palatal bars is manipulated (torque, toe-in) to avoid distortion.

## NiTi Cinch Back Pliers (060-037)



Tips are at 90° angles from the plier handle to permit access to the arch wire in making adjustments in the labial lingual direction.

## Intra Oral Finishing Pliers

(060-045,050,051,052,053)



Produce accurate first order steps (in, out, up, down) in finishing archwires without removing the archwire using this set of 4 pliers. 1/4, 1/2, 3/4 and 1 millimeter steps are created with a squeeze of the plier.



These pliers should **NOT** be used on archwires of greater size than .018.

(060-045) - Set of 4 Pliers

(060-050) - 1/4mm Pliers

(060-051) - 1/2mm Pliers

(060-052) - 3/4mm Pliers

(060-053) - 1mm Pliers

## Hole Punch Cutter (060-055)



These pliers are used on Clear aligners to punch a 6mm hole in aligner.

## Vertical Dimple Aligner Pliers - 3mm (060-057)



These pliers are used for adjusting clear aligners. Used to make a 3mm dimple anywhere on the aligner. **Not available in PDS Spain.**

## Second Order Bend Finishing Pliers (060-060,061,062,063,064)



- (060-060) - Set of 4 Pliers
- (060-061) - UR-LL-M-UL-LR-D
- (060-062) - UR-LL-D-UL-LR-M
- (060-063) - UR-LL-DZ-UL-LR-MZ
- (060-064) - UR-LL-MZ-UL-LR-DZ

These specialized pliers produce accurate second order bends (for distal or mesial root tip correction) in finishing archwires without removing the archwire from the mouth and can be used with any size wire. They produce 24 custom bends intraorally.

## Hollow Chop Pliers (060-075)



Changing the archform in rectangular archwires. These pliers are also used to curve key-hole and T-loops towards the gingival tissue to avoid lip irritation.

## Molar Bracket Holder (060-080)



This molar bracket holder has been modified to conveniently hold the first and second molar IP brackets from the buccal surface (occlusal-lingual grab) without bending the hooks.

## NiTi Chilling Instrument (060-085)

Used with 018x025 Heat Activated Wire.



## Primo Distal-End Cutter (Long Handle) (060-096)



Distal end cutter that “shears” the wire instead of “cutting” the archwire.

## IP Bracket Holder Positioners (060-100,101,102,103,104,105,110,118)



- (060-118) - Bracket Delivery Set (6 Holders)
- (060-100) - Bracket Holder - 2.5mm
- (060-101) - Bracket Holder - 3.0mm
- (060-102) - Bracket Holder - 3.5mm
- (060-103) - Bracket Holder - 4.0mm
- (060-104) - Bracket Holder - 4.5mm
- (060-105) - Bracket Holder - 5.0mm
- (060-110) - Height Verification Guide
- (090-103) - Bracket Tray Only




The Height Verification Gauge (060-110) can be placed between the two parts of the height gauge to be certain the instrument is calibrated correctly. This becomes necessary after continuous use in the practice, where the ends can be bent.

The instruments come in heights of 2.5mm to 5mm in 0.5mm increments. A minimum of 6 instruments are needed to have one instrument for each height.

Heights are dependent on the treatment selected, the curve of spee, and the esthetic setup.






## HEADGEAR

### Cervical Facebows



Produce accurate first order steps (in, out, up, down) in finishing archwires without removing the archwire using this set of 4 pliers. 1/4, 1/2, 3/4 and 1 millimeter steps are created with a squeeze of the plier.

- (070-0080) - Cervical Facebow Sz-1, 83mm
- (070-009) - Cervical Facebow Sz-2, 90mm
- (070-010) - Cervical Facebow Sz-3, 97mm
- (070-011) - Cervical Facebow Sz-4, 104mm
- (070-012) - Cervical Facebow Sz-5, 111mm



### Flexible Millimeter Ruler (019-046)

Flexible millimeter ruler used to measure the size of facebow and looped arch wires.

### Cervical Pull Neckstrap

with Force Modules


No adjustment of force is needed with this set-up. The modules apply 12oz (light modules) or 24oz (heavy modules) of force at any length past the first centimeter of extension.

- (070-018) - Headgear Modules - 12oz (2pk)
- (070-019) - Headgear Modules - 24oz (2pk)
- (070-021) - Cervical Neckstrap



### High Pull Facebows

070-005,006,007,008



Fitting the inner bow is the same with the High Pull Facebow as is with the Cervical Facebow. The outer bow is now shorter, to allow minimum adjustment for proper force application. **Available in Sizes 2-5.**

Adjust the outer bow with a 3-prong facebow plier until the “string” of the force module passes through the center of resistance.


### Hi-Pull Headcap (070-025)



Headcap: For attachment of the modules and positioning on the head of the patient.

### Reverse Pull Headgear (070-020)

Reverse pull headgear is used in Class III cases with maxillary insufficiency.





## ELASTICS



### Elastic Attachers (080-004)

Elastic Attachers are handy for patients whom lack dexterity or where placement is difficult.

Sold in packs of 100.



### 6oz 1/4" Latex-Rainbow Elastics (080-015)

6oz 1/4" elastics are colored elastics used for intramaxillary or intermaxillary applications. Sold in packs of 50.

### (080-025) - Sea Lions (6.5oz 1/4") - Latex - 50pk

Used for intramaxillary or intermaxillary applications.

### (080-030) - Eagle (4.5oz 1/4") - Latex - 50pk

Used in intramaxillary or intermaxillary mechanics. This is the elastic of choice for Class III cases.

### (080-035) - Gorilla (4.5oz 3/16", 1") - Latex - 50pk

Used for vertical applications to close the bite or finish the occlusion.

### (080-040) - Rhino (6oz 3/8") - Latex - 50pk

Used with reverse pull headgear.



### (080-045) - Marlin (6oz, 1/4") - Latex Free - 30pk

For those patients that have a sensitivity to Latex.

### (080-050) - Jellyfish (4.5oz, 1/4") - Latex Free - 30pk

Used for vertical applications to close the bite or finish the occlusion.

### (110-180) - Giant Panda (4.5oz, 5/16") - 50pk

These elastics are recommended for aligner cases that require class II elastics.

## PERIPHERAL ALIGNERS

### Chewies (110-100,101,102)

Help patients improve the fit of their aligners and close gaps, both inside and outside of the doctor's office. Sold in bags of 10 (2ea).



(110-100) - White (Unflavored)

(110-101) - Green (Mint)

(110-102) - Pink (Bubble Gum)



### Aligner Round Bondable Metal Buttons (110-105)

Sold in a package of 10.



## Aligner Round (110-110,111) Bondable Composite Buttons

Sold in a package of 5.

(110-110) - Clear Buttons

(110-110) - White Buttons

## IPR Gauges (110-115)

Verify the amount of interproximal reduction performed during a procedure with these IPR gauges.

This set contains 6 different gauges: 0.10mm, 0.20mm, 0.25mm, 0.30mm, 0.40mm, and 0.50mm.

## Combined IPR Strip Kit (110-135)

Contains 8 strips of each: yellow, red, dark blue, green; and 2 of each: clear, cyan, purple, brown.

Total of 40 strips per kit.

## Individual IPR Strips

(110-120) - Yellow (1 Side) - 0.07mm - 8pk

(110-121) - Red (2 Sides) - 0.12mm - 8pk

(110-122) - Dk Blue (2 Sides) - 0.15mm - 8pk

(110-123) - Green (2 Sides) - 0.20mm - 8pk

(110-126) - Clear (1 Side) - 0.10mm - 8pk

(110-127) - Cyan (1 Side) - 0.15mm - 8pk

(110-128) - Purple (2 Sides) - 0.25mm - 8pk

(110-129) - Brown (2 Sides) - 0.30mm - 8pk



## Perfect-A-Smile Kit (110-145,146,147,148)

Tooth colored light cure paint, the simplest way to create a tooth color pontic on any aligner, where there are missing teeth or on extraction cases. Kit includes all 3 shades.

Individual Shades Available:

(110-146) - Dark Shade

(110-147) - Medium Shade

(110-148) - Light Shade



## Outie Aligner Remover Tool (110-135)

The OUTIE Tool helps to remove aligners easily and quickly even with attachments on the teeth. Sold in packs of 5.

## WAXES



## Wax Bite Wafers with Foil Inserts (012-011)

Wax preformed in an archform for quick, accurate bite registration. Sold in packs of 24.



## Patient Wax Packs (014-001)

Essential for patients in the first month with fixed appliances, and for some, throughout treatment. Sold in packs of 50.



## PATIENT MANAGEMENT & RECORDS

### Impression Trays (012-006)

These plastic impression trays with extended borders are available in four sizes, upper and lower. Sold in packs of 25.

Sizes Available:

**Small** - Blue      **Large** - Yellow  
**Medium** - Green   **Pedo** - Red



### Model Storage Box - Black (012-010)

Each drawer stores four models. **Not available in PDS Spain.**



### Retainer Cases (013-006)

Used to help your patient protect their retainer. Sold in box of 25 assorted colors.



### Patient Hygiene Kit (014-005)

Comes with everything a patient needs for fixed appliance maintenance.



### Gray Diagnosis Form (015-001)

Organizes your diagnosis and assures notation of problems and clinical impressions seen prior to treatment. Sold individually.




### Tracing Acetate (016-002)

Used for tracing of X-rays. Sold in packs of 100.



### Protractor Template (016-012)

Used for tracing Cephalometric X-rays. **Not available in PDS Spain.**



### Truain Retainer Plastic (016-025,030)

Used for forming clear plastic retainers. Sold in packs of 25.

(016-025) - .30 Square

(016-030) - .30 Round





## Intraoral Mirror Set (100-004)

These non-fog rhodium-coated mirrors are essential in getting good pictures and thus creating the necessary records for your orthodontic cases. Set of 3. Also available individually.

(100-005) - Adult

(100-006) - Pedo

(100-007) - Side View



## Occlusal Photo Retractor (100-010)

Used to retract the lips and cheeks for the important mandibular and maxillary occlusal photographs.



## Lateral Cheek Retractors (100-015)


Used to retract the lips and cheeks for the important mandibular and maxillary occlusal photographs. Sold in pack of 2.

## PACKAGING TEMPLATES



## Metal and Ceramic Bracket WHEELS ONLY (012-023)

Empty bracket wheels with select colored label for metal and ceramic brackets are used to keep brackets organized.



## Archform Templates


(014-040,045,050)

Archform templates are clear templates used to assist you in determining the correct archwire shape. Sold individually.

(014-040) - Upper Archform

(014-045) - Lower Archform

(014-050) - Inner Archform



## Bracket Organizer (090-090)

A printed white card for all quadrants with tooth numbers to allow for single patient bracket setup



## Band Box (090-094)

Band Box with divided insert for our 34 sizes of bands. Labeled by quadrant.



## Archwire Binder Page (300-010)

Individual page inserts to keep your inventory of Upper, Lower, T-loop and Keyhole archwires organized. **Pages do not come with wires.**



## BRACKET WHEELS

**BETTER THAN ROTH INVENTORY SYSTEM**  
 with moderate to severe lower crowding  
 when severe lower crowding is present  
 La/Sla  
 Ne

### Archwire Binder with Page Inserts (090-097,098,099,100)



Binders with page inserts with to keep your inventory of archwires. Each insert page corresponds to one size of archwire and can keep all different forms. **Pages do not come with wires.**

- (090-097) - Upper Archwire
- (090-098) - Lower Archwire
- (090-099) - T-Loop Archwire
- (090-100) - Keyhole Archwire



### Patient Box (090-106)

Patient boxes are a great way to stay organized and store individual patient appliances.



### Ortho Model Holders (014-060)

Used to position the upper and lower study models in the closed bite position when scanning. Sold in sets of 5. **Not available in PDS Spain.**

### Complete Metal Bracket Repair Kit (011-004)

The Complete Metal Bracket Repair Kit consists of 10 bracket wheels (2 of each) to be used as back-ups in the clinic in case a patient loses brackets, or you need to alter your prescription. Individual wheels sold in packs of 2.



#### Individual Wheels Available:

- (011-0030) - Peach - Super Lingual Torque
- (011-0031) - Brown - Super Labial Torque
- (011-0041) - Orange - Incisor Labial Torque
- (011-0042) - Burgundy - Incisor/Cuspid Rotation
- (011-0043) - Yellow - Bicuspid Rotation
- (011-0044) - Violet - Incisor/Cuspid Lingual Torque
- (011-0045) - Light Purple - Standard Variations
- (011-0046) - Pink - Molar Variations
- (011-0047) - Blue - Molar/Canin/Incisor Variations
- (011-0048) - Green - Torque Rotation Variations

### AM Metal Bracket Wheel Kit (011-0060)

Absolute Minimum Kit comes with all 10 wheels (1 of each).



#### AM Bracket Wheels sold individually:

- (011-0050) - Peach - Super Lingual Torque
- (011-0051) - Brown - Super Labial Torque
- (011-0052) - Orange - Incisor Labial Torque
- (011-0053) - Burgundy - Incisor/Cuspid Rotation
- (011-0054) - Yellow - Bicuspid Rotation
- (011-0055) - Violet - Incisor/Cuspid Lingual Torque
- (011-0056) - Light Purple - Standard Variations
- (011-0057) - Pink - Molar Variations
- (011-0058) - Blue - Molar/Canin/Incisor Variations
- (011-0059) - Green - Torque Rotation Variations



## Ceramic Bracket Wheels (011-006,007)

State-of-the-art ceramic brackets are strong and durable with nearly an invisible aesthetic. Sold as a convenient wheel.

**(011-006)** - Green - La/Li Variations

**(011-007)** - Purple - Standard Variations



## Standard Bracket Set - U/L 5-5 (019-011)

Standard Roth Bracket Set Upper and Lower 5-5, one of each.

## Better Than Roth Kits

The Better Than Roth (BTR) system allows you to follow the step-by-step POS treatment plans while having the perk of our bulk discount to build inventory for your practice.



### BTR Metal Color Kits Available:

**(018-002)** - Green - Variable Torque Class II

**(018-006)** - Blue - Variable Torque Class 1 NE

**(018-011)** - Red - Variable Torque Class II NE Class II Elast.

**(018-016)** - Yellow - Variable Torque Class I or II Bi Ext.

**(018-017)** - Orange - Class II Div II Non-Extraction

**(018-018)** - Gray - Class III Non-Extraction

### BTR Ceramic Color Kits Available:

**(018-020)** - Aqua - Ceramic Class II, Non-Ext. (4mm or less)

**(018-021)** - Pink - Ceramic Non-Ext. Distilization Up. Molar

**(018-022)** - White - Ceramic Class I, II, III Mod. Crowding

**(018-023)** - Purple - Ceramic Bicuspid Ext. Class I & II

Purchase any 50 BTR Metal brackets for 35% off or 10 BTR Ceramic brackets for a 10% off discount! Use code **BTR METAL** or **BTR CERAMIC** at checkout.  
(Cannot mix Metal & Ceramic)

## BAND FITTING KITS

### Complete Band Fitting Kit (015-010-00)

The Complete Band Fitting Kit comes with molar fitting bands for all 4 quadrants, along with the upper and lower bicuspid. The 4 molar quadrants contain 2 fitting bands of every size from 1-34 and 1 band each of upper and lower bicuspid sizes from 1-24.



#### Individual Band Fitting Kits (1/ea):

**(015-014)** - Upper Bicuspid Band

**(015-016)** - Lower Bicuspid Band

**(015-017)** - UR 1st & 2nd Molar

**(015-018)** - UL 1st & 2nd Molar

**(015-019)** - LR 1st & 2nd Molar

**(015-020)** - LL 1st & 2nd Molar

### AM Band Fitting Kit (015-010-001)

Absolute Minimum Band Fitting Kit comes with all 4 quadrants containing 1 molar fitting band of every size from 1-34.



#### Individual Band Fitting Kits (1/ea):

**(015-010-002)** - UR 1st & 2nd Molar

**(015-010-003)** - UL 1st & 2nd Molar

**(015-010-004)** - LR 1st & 2nd Molar

**(015-010-005)** - LL 1st & 2nd Molar

## SKELETAL ANCHORAGE



### 17x25 Stainless Steel Archwire (030-0065)

A tube of 17x25ss archwire is needed to make sectional wires with helices for securing skeletal anchorage to the IP Appliance "pig" molar buccal tube variation. Sold in tubes of 10.




## 6mm NiTi Closed Coil Springs (040-019A)



These coils were specially made for skeletal anchorage applications where there is limited distance to activate a coil.

## Skeletal Anchorage Kit (120-005)

This pre-assembled kit contains the autoclavable container with a small inventory of bone plates, screwdriver, driver tip with screw securing sleeve, contra-angle driver tip, and ceramic tray for storing bone screws and ortho anchors.



## Bone Truscrews (120-069,070,071)

Bone screws are used to secure bone plates. Sold in packs of 5. **Not available outside USA.**

(120-069) - 1.2 x 4.5mm

(120-070) - 1.5 x 6mm

(120-071) - 1.5 x 7.5mm



## Ortho Anchors (120-075,076)

Ortho anchors are simply a bone screw with a larger head to allow securing of ligature ties and coils. Sold individually. **Not available outside USA.**

(120-075) - 1.5 x 6mm with cuff

(120-076) - 1.5 x 8mm with cuff



## Tissue Punches - 2.5mm (120-085)

When placing an orthodontic anchor into attached gingiva on the buccal or palate, a soft tissue (biopsy) punch is first used to gain access to the bone before drilling a pilot hole. Disposable. Sold in packs of 10. **Not available outside USA.**

## Spool of .0126 SS Wire 28G (120-090)

It is recommended that .012 ligature ties be used to secure coils to the bone plates and archwires. .012 ligature is provided on a spool that may be autoclaved. **Not available outside USA.**



## Gram Gauge (120-095)

The gram gauge is used for measuring the force of the coils used in Skeletal Anchorage. (25-250 grams)



## Bone Plates (120-100,101,102,103,080)

The gram gauge is used for measuring the force of the coils used in Skeletal Anchorage. Sold individually. **Not available outside USA.**


(120-080) - T-Shape - 7 hole (18mm)

(120-100) - Y-Shape - 6 hole (31mm)

(120-101) - 8 hole (31mm)

(120-102) - 4 hole (15mm)

(120-103) - 2 hole (7mm)





### Pilot Twist Drill (120-105,106)

Straight drills of 1.1mm and 0.8mm are for drilling a pilot hole before attempting to place the bone screw or ortho anchor. **Not available outside USA.**

- (120-105) - 1.1mm
- (120-106) - 0.8mm

### 1.25mm Driver Tip (120-107,108)

The driver tip is placed into the screwdriver handle after pulling back on the end of the handle and releasing to lock. The driver tip can therefore be replaced without the need for replacing the handle. **Not available outside USA.**

- (120-075) - 55mm
- (120-076) - 19mm

### Bone Screw Driver Handle (120-110)

Specific Driver handle for our self-drilling Bone Screws. Sold individually. **Not available outside USA.**

### Jeil Micro Auto Screws (130-125,126)

Made with Titanium Alloy. Available in two sizes. Sold individually.

- (130-125) - 5mm
- (130-126) - 6mm

### Jeil Micro Straight Plate - 8 Hole

(130-130)

Micro plate thickness 0.5mm. Hole diameter 01.5mm. Made with titanium.

### Jeil Mid Double Y Plate - 6 Hole

(130-135)

Plate thickness 0.6mm. Width 8.9mm. Made with titanium.

### Jeil Screw Driver Handle - Long

(130-140)

Screwdriver handle is used to drive screws into the bone. Made with Aluminum.

### Jeil Micro Screw Driver Shaft

(130-145)

Length 70mm. Used with the Jeil Screwdriver Handle. Made with stainless steel.

### Jeil Dual Top Anchor TAD Hexed Driver Shaft

(130-155,156)

Driver shaft is only compatible with the Jeil Dual Top Anchor TAD Bone Screws (6mm or 8mm).

- (130-155) - 20mm
- (130-156) - 70mm



## Jeil Dual Top Anchor TAD Bone Screws (130-165,166)

Length 1.4mm. Self-tapping, immediate loading, easy removal. Made with titanium alloy. Sold individually.

(130-165) - 6mm

(130-166) - 8mm

## EDUCATION & OFFICE KITS



### IP Seminar Kit (019-002)

The IP Seminar Kit is a group of products that will be used with all hands-on exercises (such as cements, brackets and ligature wires). **This kit cannot be modified and does not include any pliers.**



### Educational Kit (Kit-E)

The Educational Kit has all the products you will need to participate in the hands-on portion of class.

The Ed-Kit consists of two parts: the IP Seminar Kit and any additional instruments needed for in-class exercises.

## Absolute Minimum Kit (Kit-AM)



Absolute Minimum Kit is designed to help you choose the product that is absolutely essential to you being able to create an orthodontic department in your practice. Eventually, you will need all of the products contained in the Practice Kit, which would be a better step in the right direction to completely set up your practice more economically.



## Practice Kit (Kit-P)

The Practice Kit was designed to help you set your practice up correctly to easily incorporate orthodontics. This kit is optional, however we believe that this kit allows you to properly introduce orthodontics quickly and easily, and encourage you to purchase this to set up your practice for maximum efficiency.




## Aligner+ Kit (Kit-A+)

The Aligner+ Kit contains the products needed to participate in the hands on exercises for Seminar 2 of the Aligners+ Series. These products will also be used in the practice.



## BOOKS, BROCHURES & OFFICE MARKETING




### Braces Owner Manual (014-006)

A brief guide to the wearing and care of braces, with illustrations and brief descriptions of the common uses of orthodontics. Also sold as part of the patient hygiene kit.



### Elastics Brochure (014-007)


Brief explanation of why and how to wear elastics successfully.



### Orthodontic Brochures (014-018,021,022)

A guide for adults, teens and kids explaining Orthodontics, why it's necessary and how to care for their appliances.

- (014-018) - Adult Brochure
- (014-021) - Teen Brochure
- (014-022) - Early Orthodontics Treatment



### Contemporary Orthodontics Author: William Proffit (014-023)

Good book for the student who is knowledgeable in orthodontics. **Not available in PDS Spain.**



### POS Patient Brochures (014-024,025,026)


Explains the Individual Patient (IP) Appliance and why it is the most effective and efficient way to straighten teeth.

- (014-024) - Teen Brochure
- (014-025) - Ortho for Children Brochure
- (014-026) - Teen Brochure (Spanish)



### Tabletop Brochure Holder (014-027)

Place your "POS Patient Brochures" in a brochure holder to attract people's attention in your lobby or consultation room.



### DA Illustrated Handbook for Ortho Assisting (014-205)

This handbook is meant to serve as a manual that can be tailored to meet the unique needs of any office. Comes with customizable forms to allow you to personalize each section according to your individual practice.

## Orthodontic Office Posters (014-085,200,201)

Orthodontic office posters are a great way to let your patients know that you now offer orthodontics. **Not available in PDS Spain.**



**A great smile is part of your style.**  
Our office uses the Individual Patient (IP) Appliance® for faster results.  
ASK YOUR DENTIST  
For the most advanced braces from the United States can improve your smile!



**GET THE SMILE OF YOUR DREAMS!**  
with Advanced, Computer-Customized Braces  
ASK YOUR DENTIST  
For the most advanced braces from the United States can improve your smile!



**Give it your all and get ready for a winning smile.**  
Our office uses the Individual Patient (IP) Appliance® for faster results.  
Orthodontics fits your kid's active lifestyle, ask us how.

- (014-085) - Teen Office Poster
- (014-200) - Young Adult Office Poster
- (014-201) - Young Kids Office Poster

## Social Media Marketing Package (014-204)

Be successful at marketing your orthodontic services to your patients by maximizing your social media opportunities with these pre-designed posts. **Not available in PDS Spain.**

- Digital Assets:
- 6 Social Media Posts (plus 2 bonus posts, one for Pinterest and LinkedIn)
  - 1 Content Planning Calendar
  - 1 Website (individualortho.com)



Note: Social media posts are different for the General Practice Marketing Package and the Social Media Marketing Package



## General Practice Marketing Package (200-200)

The General Practice Marketing Package contains materials that will help you successfully market your orthodontic services to your patients to build new business. **Not available in PDS Spain.**

Physical components:

- 2 Posters - 1 Teen / 1 Kid
- 20 Patient Brochures - 10 Teen / 10 Kid
- 2 Brochure Holders
- 1 USB Drive with Monitor Slides
- 3 Consultative Treatment Guides



Digital Assets:

- 10+ Social Media Posts & Instructions
- 1 Content Planning Calendar
- 1 Website (individualortho.com)



**ALL THIS IN ONE GREAT PACKAGE!**



**Give it your all and get ready for a winning smile.**  
Our office uses the Individual Patient (IP) Appliance® for faster results.  
Orthodontics fits your kid's active lifestyle, ask us how.

TYPES OF ORTHODONTIC CASES			
<b>CHROMIUM BETH</b>	<b>WIDELY SPACED TEETH</b>	<b>OVERBITE</b>	<b>UNDERBITE</b>
<ul style="list-style-type: none"> <li>• Improved oral health</li> <li>• Better appearance</li> <li>• Better function</li> </ul>	<ul style="list-style-type: none"> <li>• Close gaps between teeth</li> <li>• Better oral hygiene</li> <li>• Better function</li> </ul>	<ul style="list-style-type: none"> <li>• Reduce jaw strain</li> <li>• Improved oral hygiene</li> <li>• Better function</li> </ul>	<ul style="list-style-type: none"> <li>• Better oral function</li> <li>• Better oral hygiene</li> <li>• Better function</li> </ul>
<b>OPEN BITE</b>	<b>CLOSED BITE</b>	<b>CROSS BITE</b>	<b>Ask your dentist if orthodontics is right for you!</b>
<ul style="list-style-type: none"> <li>• Better oral function</li> <li>• Better oral hygiene</li> <li>• Better function</li> </ul>	<ul style="list-style-type: none"> <li>• Reduce jaw strain</li> <li>• Better oral hygiene</li> <li>• Better function</li> </ul>	<ul style="list-style-type: none"> <li>• Reduce jaw strain</li> <li>• Better oral hygiene</li> <li>• Better function</li> </ul>	



**A great smile is part of your style.**  
Our office uses the Individual Patient (IP) Appliance® for faster results.  
Ask us how orthodontics will raise your selfie WOW factor.

# PROGRESSIVE DENTAL SUPPLY

135 Columbia, Suite 101  
Aliso Viejo, CA 92656  
US: +1-714-973-2266  
EU: +34-94-6361097



*Scan QR Code To Shop Now At*  
***www.SmileStream.com***